

GATE

REHAB DEVELOPMENT AB™



Hygiene chairs

Accessories

Accessories

Adjustable backrest/foam
56-258 (BASIC/DOUBLE/XL)
 For improved comfort/hygiene. The backrest provides the patient with greater depth. Machine washable at 160° F.

Safety bar
56-212- (STD/XL)
 For safer sitting. Prevents the patient from sliding out of the chair.

Side supports
56-216 For comfortable sitting. Provides soft support for the smaller individual.

Soft seat with hygiene cut-out
 (BASIC/DOUBLE)
56-206 Increased patient comfort.

Plug for soft seat
56-232 Used together with **56-206** it provides a soft shower seat for more comfortable showers.

Widening kit
56-221 Widens armrests by 2" per pair. i.e. 1" per side. For patients who need extra space.

Adjustable angle amputation support
 AMFIBI
56-228-L left
56-228-R right
 BASIC
56-289-L left
56-289-R right
 Adjustable by length and angle.

Guide
56-218 For installing waste container. Also fits standard

Adjustable angle calf support
 AMFIBI
56-227-L left
56-227-R right
 BASIC
56-288-L left
56-288-R right
 Adjustable by length and angle.

Waste container with lid
56-209 Large handle for easy handling.

Directionally locking castor
 (DOUBLE)
56-359 49"
56-355 39"
57-043 29"

Castor with brake
 (DOUBLE)
56-251 49"
56-222 39"
56-387 29"

56-257-L left
56-257-R right
Footrest for adjustable-angle leg support
 (DOUBLE/XL)
 Used in combination with amputation stump support or calf support.



Drive wheel kit
56-215 (22")
 Only fits BASIC. Incl. brakes with quick-release function so that patients can move around in hygiene chairs by themselves.



Splash guard
56-213
 The guard prevents urine leaks at the front.



Calf support
56-225
 Provides extra support for the lower leg.



Shower seat
56-205
 (BASIC/DOUBLE)
 Seat with small hole for increased comfort when showering.



TR EQUIPMENT

TR GROUP, Inc., 2990 N. Peryville Rd., Suite 2300, Rockford, IL 61107
 Phone: 800 752 6900 Fax: 815 633 0549

Hygiene Equipment Designed for People

www.trequipment.com