Height Adjustable Battery **Operated Shower Chair**



Generic Tender Specifications for Height Adjustable Battery Operated Shower Chair

- Equipment shall be delivered 100 % assembled and ready for use after charging the battery •
- Equipment is electrically height-adjustable and tiltable for caregiver's ergonomic use and is specifically • designed and manufactured for safe in-patient hospital use
- Tilting range from 25 degrees reclining to 7 degrees forward
- Shower chair chassis design to create possibilities for toileting over standard toilets •
- Handheld control unit for height settings, recline and forward tilting, for caregiver convenience •
- Shower chair to be designed for smooth and safe patient transfers from sling hoists and • standing/raising aids
- Shower chair to include a hip belt for patient safety •
- Shower chair to be equipped with an adjustable headrest for head and neck comfort, support and security ٠
- Foldable armrests with grip knobs for residents comfort and safety •
- Shower chair equipped with emergency stop
- Shower chair equipped with four twin castors with brakes for safe and easy patient transport
- 24 V removable and rechargeable batteries (min. 5 Ah) to be charged via main supply and with visible low battery indicator
- All electronics on the Shower Chair to be in the protection class IPX6 or higher •
- Minimum of 40 up-down movements with fully charged batteries with a 150 kg/330 lbs. load •
- Seating surfaces shall be made of a soft, durable material that is non water absorbent, easy to clean, disinfect and comfortable to sit on
- Anti-corrosion treatment on all sheet metal parts

Specification text to add for options/accessories if needed

- Shower chair to include a bedpan (lid included)
- Shower chair to be equipped with adjustable and foldable footrest for security and comfort for residents •
- Shower chair to include front closing armrest for extra patient safety •
- Shower chair to include a chest belt for additional security and support for patient

Dimension and load specifications

- Weight Limit (patient): 200 kg/440 lbs. Operating Height Range: • Width Seating Area: • Depth Seating Area:
- External Width:
- **Castor Diameter**

Phone: +46 140 38 50 60

info@trequipment.se www.trequipment.com Min: 535 mm/21 inch Min: 500 mm/19^{2/3} inch Min: 430 mm/16^{4/5} inch Max: 700 mm/271/2 inch Min: 100 mm

Max: 935 mm/364/5 inch Max: 550 mm/21^{2/3} inch Max: 450 mm/17^{2/3} inch



Group Subsidiaries TR Equipment AB, Sweden

TR Group Inc, US, Phone: +1 (800) 752 6900 TR Equipment Ltd, UK, Phone: +44 (0) 844 335 8386 TRanspoRehaMed GmbH, Germany, Phone: +49 (0) 212 645850-0



