

GATE

REHAB DEVELOPMENT AB™



Hygiene chairs

Amfibi XL

Amfibi XL

EL

Amfibi XL was designed for the bigger, heavier patient and therefore provides extra space and strength and electrical height adjustment. Maximum user weight is 550 pounds.

Can handle
550 lbs!

Plenty of room

Maximum user weight for the Amfibi XL is 550 lbs and seat depth is larger than standard. The hygiene chair has also been widened to 26" between armrests to provide extra room and good stability.

Durable materials

Manufactured from corrosion-proof materials.

Lockable armrests

Amfibi XL has armrests that can be flipped up and locked in position to make it easier for the patient to be of help with the transfer.

Easy control

Electrical adjustment via a cabled remote which is sealed to a high protection class (IP65) for optimal safety.

Built-in soft seat

The XL's seat and back are adapted for an upright posture that does not close the pelvis. The soft, ergonomic seat (PU) is standard and built in to the chair. It minimizes pressure on the hips, tendons and blood vessels and is designed to suit both men and women.

Adjustable back-strap

Adjustable backrest The soft, comfortable back-strap relieves pressure on bony areas. It is easily adjustable for best support and seat depth. Machine washable.

Easy to clean

The seat is easy to remove and clean, while its rounded front edge contributes to a comfortable leg position.

Directionally locking/lockable castors

One of the castors can be locked directionally for improved safety during transit. The other castors are lockable.

Swing-away, detachable footrests

The footrests flip-up, swing-away and are easily detachable.

Electric height adjustment

Amfibi XL has a unique pillar design for quick and easy seat height adjustment. The seat can be raised by up to 32" from its lowest position and thus requires minimum effort from care personnel.

Guides and waste containers are available as accessories.

Technical specifications

CE CE marked according to MDD 93/42/EEC.

Part No.	Product	Width, external	Max user weight	Seat height	Width between armrests	Seat depth
56-204	Amfibi XL	31"	550 lbs	22" - 32"	26"	19"



TR EQUIPMENT

TR GROUP, Inc., 130 North Madison Street, Rockford, IL 61107
Phone: 800 752 6900 Fax: 815 633 0549

Hygiene Equipment Designed for People

www.trequipment.com