



TR 1700

Height adjustable bathtub

Hygiene Equipment
Designed for People

Quick guide



1	Thermostatic mixing valve
2	Shower handle
3	Overflow and drainage locking
4	Aluminium actuator column
5	Adjustable feet
6	Hand control for height adjustment
7	Actuator

For further information see Operating manual

Working Load Limit: 150 kg (330 lbs)

Head office

TR Equipment AB, Sweden
Phone: +46 140 38 50 60
info@trequipment.se
www.trequipment.com

Group Subsidiaries

TR Group Inc, US, Phone: +1 (800) 752 6900
TR Equipment Ltd, UK, Phone: +44 (0) 844 335 8386
TRanspoRehaMed GmbH, Germany, Phone: +49 (0) 212 645850-0



TR EQUIPMENT