

# Bonn

## Shower and Commode Chair

*Hygiene Equipment  
Designed for People*

### Rebotec Bonn

A durable, long lasting commode and shower chair with an easy-to-clean design.



Registered at Medical Products  
Agency, Germany,  
as a MD Class 1 product,  
MDR 2017/745.

TR reserves the right to make  
design changes and/or  
improvements without notice.



Rebotec Bonn has soft and durable polyurethane seat cover for comfortable seating and easy cleaning. The shower chair has 5" hospital grade swiveling and locking castors for secure transfer. Hygiene opening can easily be moved from the front to the back. Hygiene bucket with vacuum seal top included.



**TR EQUIPMENT**





# Bonn

## Shower and Commode Chair

Hygiene Equipment  
Designed for People



### Dimensions



### Technical specifications

Dimensions		
A. Height	101 cm	40"
B. Width	56 cm	22"
C. Depth	87 cm	34"
D. Seat height	53 cm	20,8"
E. Seat width	46 cm	18"
F. Seat depth	45 cm	18"
G1. Max. height toilet (with pan receptacle)	43 cm	17"
G2. Max. height toilet (without pan receptacle)	47 cm	18,5"
H. Foot height	33-48 cm	13"-19"
I. Castor	13 cm	5"
Body weight max.	130 kg	290 lbs
Product weight	11,9 kg	26,2 lbs

### Equipment

- Ergonomically shaped backrest with integrated hand grip for caregiver aid
- Swinging armrests enable easy side transfer
- Height adjustable, swinging and removable foot supports standard
- 5" hospital grade swiveling and locking castors for secure transfers
- Durable polyurethane seat cover allows for easy cleaning
- Polyamide nylon construction guarantees no rust
- Hygiene bucket with vacuum seal top included
- Approved up to 290 lbs body weight
- Five year warranty on frame for long lasting confidence

### Accessories

- Safety Belt for waist
- Belt Comfort Pad
- Arm rest
- IV Pole

### Art no

322.30.20  
322.10.50  
440.10.32  
440.00.93

#### Head office

TR Equipment AB, Sweden  
Phone: +46 140 38 50 60  
info@trequipment.se  
www.trequipment.com

#### TR Equipment UK Ltd

Innovation Centre, Gallows Hill,  
Warwick, CV34 6UW. England  
Phone: +44 (0) 844 335 8386  
Email: info@trequipment.co.uk  
Site: www.trequipment.co.uk



TR EQUIPMENT