



Bonn

Shower and Commode Chair

Hygiene Equipment
Designed for People

Rebotec Bonn

A durable, long lasting commode and shower chair with an easy-to-clean design.



Registered at Medical Products Agency, Germany, as a MD Class 1 product, MDR 2017/745.

TR reserves the right to make design changes and/or improvements without notice.



Rebotec Bonn has soft and durable polyurethane seat cover for comfortable seating and easy cleaning. The shower chair has 5" hospital grade swiveling and locking castors for secure transfer. Hygiene opening can easily be moved from the front to the back. Hygiene bucket with vacuum seal top included.



TR EQUIPMENT





Bonn

Shower and Commode Chair

Hygiene Equipment
Designed for People

Dimensions



Technical specifications

Dimensions		
A. Height	101 cm	40"
B. Width	56 cm	22"
C. Depth	87 cm	34"
D. Seat height	53 cm	20,8"
E. Seat width	46 cm	18"
F. Seat depth	45 cm	18"
G1. Max. height toilet (with pan receptacle)	43 cm	17"
G2. Max. height toilet (without pan receptacle)	47 cm	18,5"
H. Foot height	33-48 cm	13"-19"
I. Castor	13 cm	5"
Body weight max.	130 kg	290 lbs
Product weight	11,9 kg	26,2 lbs

Equipment

- Ergonomically shaped backrest with integrated hand grip for caregiver aid
- Swinging armrests enable easy side transfer
- Height adjustable, swinging and removable foot supports standard
- 5" hospital grade swiveling and locking castors for secure transfers
- Durable polyurethane seat cover allows for easy cleaning
- Polyamide nylon construction guarantees no rust
- Hygiene bucket with vacuum seal top included
- Approved up to 290 lbs body weight
- Five year warranty on frame for long lasting confidence

Accessories

- Safety Belt for waist
- Belt Comfort Pad
- Arm rest
- IV Pole

Art no

322.30.20
322.10.50
440.10.32
440.00.93

Head office
TR Equipment AB, Sweden
Phone: +46 140 38 50 60
info@treqquipment.se
www.treqquipment.com

TR Group Inc
130 North Madison Street
Rockford, IL61107, USA
Phone: +1 (800) 752 6900
Email: information@treqquipment.com
www.treqquipment.com



TR EQUIPMENT