

Dallas Bariatric Shower and Commode Chair

Rebotec Dallas

A durable, long lasting bariatric shower chair with an easy-to-clean design. Polyamide nylon construction with stainless steel fasteners guarantees no rust.







Registered at Medical Products Agency, Germany, as a MD Class 1 product, MDR 2017/745.

TR reserves the right to make design changes and/or improvements without notice.


It has soft and durable polyurethane seat for comfortable seating and easy cleaning. Removable rear support and seat cover for toileting. Shower-commode chair with 5" castors. Hygiene opening can easily be moved from the front to the back. Hygiene bucket with vacuum seal top included.





Dallas Bariatric Shower and Commode Chair

Dimensions



Technical specifications

Dimensions

A. Height

B. Width

C. Depth

D. Seat height

E. Seat width


F. Seat depth

G1. Max. height toilet

G2. Max. height toilet

(with pan receptacle)

(without pan receptacle)



Equipment

Ergonomically shaped backrest with integrated hand grip for caregiver aid

- Swinging armrests enable easy side transfer
- Hygiene opening can easily be moved from the front to the back
- Height adjustable, swinging and removable foot supports standard
- 5" hospital grade swiveling and locking castors for secure transfers
- Soft and durable polyurethane seat allows for comfortable seating and easy cleaning
- Polyamide nylon construction with stainless steel fasteners guarantees no rust
- Hygiene bucket with vacuum seal top included
- Approved up to 440 lbs body weight
- idence.

H. Foot height	38 - 53 cm	15" - 21"	 Five year warranty on frame for long lasting confidence 	
I. Castor	13 cm	5"		
Body weight max.	200 kg	440 lbs	AcessoriesSafety Belt, for chest and waistBelt Comfort Pad	Art no 322.30.20 322.10.50
Product weight	22 kg	48 lbs		
			Arm rest	440.10.32
			 Splash-Protection 	440.10.40

Head office

TR Equipment AB, Sweden Phone: +46 140 38 50 60 info@trequipment.se www.trequipment.com

TR Group Inc

104 cm

74 cm

90 cm

56 cm

64 cm

48 cm

45 cm

48 cm

41"

29"

22"

25"

19"

17,7"

19"

35,4"

130 North Madison Street Rockford, IL61107, USA Phone: +1 (800) 752 6900 Email: information@treuipment.com www.trequipment.com

