



TR Spaceline

Height Adjustable Bathtub

Hygiene Equipment
Designed for People

TR Spaceline - Height Adjustable Bathtub

The TR Spaceline is a freestanding, height-adjustable bathtub designed to provide a relaxing and satisfying bathing experience. Its extra width accommodates plus-size users and offers additional space for users who want to practice movement within the tub. A range of sensory features is available, including an LED lighting system, a music soundwave system with four speakers (AUX and Bluetooth), and a 12-jet air-spa system.

The Spaceline's straightforward installation allows it to be placed in various central bathing areas or bathrooms. The controls are designed for ease of use, and its sleek design ensures simple cleaning and maintenance.



Registered at Medical Products Agency, Sweden as a MDR Class 1 product according to MDR 2017/745.

TR reserves the right to make design changes and/or improvements without notice.



Specifications

Length:	2116 mm / 83"
Width:	936 mm / 36.8"
Working Load Limit:	250 kg / 551 lbs
Typical lift speed, no load:	8 mm/s / 0.3"/s
Operating height range:	650 - 1000 mm / 26" - 40"
Power supply:	230 V, 50 Hz, 10 A / 110 V, 60 Hz
Overall weight (empty):	115 kg / 253 lbs





TR Spaceline

Height Adjustable Bathtub

Hygiene Equipment
Designed for People

Features

- Innovative and attractive design
- Easy to operate and learn - Quick introduction
- Easy access to the bathtub from three sides
- Ergonomic tub shape
- Easy to clean and maintain
- Detached, no wall or floor fixing

Standard Equipment

- High quality bathtub with side covers of reinforced fibreglass
- Hand shower with long hose
- Overflow control and pop-up drainage incl drainage hose
- White as standard color
- Two chrome grab handles
- Removable foot support, adjustable in three positions
- Reinforced hoses for water installation
- Two wide stainless steel legs with adjustable feet



Options/Accessories available for Spaceline

- Cleaning/disinfection system
- Autofill with 2 different pre-set levels
- Bathtub and bathtub side covers in optional RAL - Color (Specify RAL code on order)
- Reverse installation configuration kit with extended hoses
- FMM Antibactus mixer (Legionella)
- Mixer Grohe 2000, with Grohe CoolTouch system

Wellness Options

- Air Spa from twelve integrated jets stimulate blood circulation and a calming environment for relaxation
- LED Color Light System with light that changes the water into different colors
- Music sound wave system with four speakers AUX and bluetooth

Application Examples

Integrated with the TR 9650 Mobile Patient Combilift the system provides safe moving and handling for both patient and attendant. During the bathing cycle the client never leaves the comfort and security of the chair or stretcher.

Note! Working Load limit for TR 9650 is not equal to working load limit for TR Spaceline.



Art. no.	Model
05196000	Spaceline 40mm drainage 230V
05196002	Spaceline 50mm drainage 230V
05196005	Spaceline 40mm drainage TMV3 230V (UK Version)
05196041	Spaceline 400 mm drainage TMV3 230 V (UK Version)

Dimensions (mm)

